# ENTRAL<sup>TM</sup>

**SONDE NASO-GASTRIQUE DE GAVAGE**  
**NASOGASTRIC ENTERAL FEEDING TUBE**  
**SONDA D'ALIMENTAZIONE ENTERALE NASO GASTRICA**  
**SONDA NASOGÁSTRICA PARA ALIMENTACIÓN FORZADA**  
**ENTERALE NASALE MAGENSONDE**

## **Description :**

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## **Description:**

The ENTRAL Paediatric naso-gastric feeding tube are designed for continuous or intermittent enteral feeding.  
Radio-opaque medical grade polyurethane tube or with radio opaque stripe.  
The graduations in centimetres printed on the tube to make insertion and checking of migrations easier..  
TUBE EQUIPPED WITH A SAFETY CONNECTOR LUER LOCK MALE in line with EN 1615 European Standard.

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## POSITIONING INSTRUCTIONS

Note: The following instructions are only recommendations for the use of this device. Specific situations may require amendments to this procedure. Determine the tube length needed to enter the stomach. Use the markings on the tube to avoid excessive length. Insert the tube through one of the nostrils. Slide the tube through the oesophagus into the stomach. Should the patient cough during this step of the procedure, the tube may have entered the trachea. Withdraw the tube immediately and start the procedure over again.

## PRECAUTIONS BEFORE USE

Systematically check the position of the tube:  
- either by sampling gastric juice (using a 50 or 60 ml syringe)  
- or by using a stethoscope during air injection (with a 50 or 60 ml syringe).  
For high risk patients, the tube's position needs to be confirmed through radiology.  
Inserting the tube in the trachea, bronchi or lungs may entail serious risks for the patient.  
Always check that the tube does not coil in the patient's back throat.  
Never use a syringe with a less than 50 ml content when checking the tube's position or rinsing. Smaller syringes create higher pressure and may seriously damage the tube and therefore harm the patient.  
Rinse the tube with 20 or 30 ml of water (using a 50 or 60 ml syringe) at every feed interruption and systematically every 6 hours.  
Should the tube be blocked, do not attempt to unblock it. It is strongly recommended to replace a malfunctioning tube and strictly forbidden to use a rod or any other device to unblock it.  
Change the device every 3 to 4 weeks.  
Single use only. Do not re-use or re-sterilise the device.

Volume intérieur approximatif des sondes ENTRAL™ pédiatriques.  
Approximative ENTRAL™ feeding tube volume  
Volumen interior de las sondas ENTRAL™ pediátricas.  
Innengröße der pädiatrischen ENTRAL™ Sonden.

	Ch 4	Ch 6	Ch 8	Ch 10
Lg 40 cm	0.4 ml	0.5 ml	0.7 ml	1.4 ml
Lg 50 cm	0.5 ml	0.7 ml	0.9 ml	1.7 ml
Lg 80 cm	0.8 ml	1.1 ml	1.4 ml	2.8 ml
Lg 120 cm	1.1 ml	1.6 ml	2.1 ml	4.2 ml

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