ENTRALTM VAC

SONDE D'ASPIRATION NASO-GASTRIQUE A DOUBLE COURANT DE TYPE SALEM
DOUBLE LUMEN NASOGASTRIC SUCTION TUBE TYPE SALEM
SONDA DOPPIO LUME NASO GASTRICA D'ASPIRAZIONE TIPO SALEM
SONDA DE ASPIRACIÓN NASOGÁSTRICA DE DOBLE CANAL TIPO SALEM.
NASO GASTRALE 2-LUMIGE ABSAUGSONDE TYP SALEM

Description :

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Description:

The ENTRAL VAC tubes are designed for gastric suction (nasogastric or nasointestinal). The double lumen allows the ability to limit depression and to protect the stomach wall. The tubes are of biocompatible polyurethane and radiopaque or have a line which can be localised by X-rays.

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INSTRUCTIONS FOR USE

Warning: The following procedure is intended for reference only, particular situations may require modifications.

1. Explain procedure to patient.
2. Place the patient in sitting position, leaning slightly backwards.
3. Measure tube for placement. Extend tube from the patient's earlobe to the bridge of the nose and from there to the tip of the xiphoid process. Reference marks are located 25, 50 and 75 cm from the distal end of the tube.
4. ENTRAL VAC tubes are prelubricated and do not need any particular preparation before insertion. It is possible to improve lubrication by applying some silicone gel on the tube.
5. Insert the tube following hospital procedure. Swallowing makes insertion easier (with or without liquid absorption).
6. Check the position of the tube by:
- X-ray control - fiber optics - aspiration of stomach juice - stethoscope during air injection
7. Connect tube to vacuum source and start suction.

CAUTION

Single use product, do not reuse. Use only syringes with volumes of at least 50 cc's to check the correct position. Replace tube after 3 to 4 weeks.

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